

PLACA BLOQUEADA LPS ÚMERO PROXIMAL EM AÇO

*LOCKING PLATE LPS PROXIMAL
HUMERUS IN STAINLESS STEEL*

*PLACA BLOQUEADA LPS HÚMERO
PROXIMAL EN ACERO*



ANVISA 80062900021

PLACA BLOQUEADA LPS ÚMERO PROXIMAL EM AÇO

LOCKING PLATE LPS PROXIMAL HUMERUS IN STAINLESS STEEL

PLACA BLOQUEADA LPS HÚMERO PROXIMAL EN ACERO



Tem como indicação, finalidade e/ou uso a estabilização e fixação de diversos tipos de fraturas, proporcionando a recuperação de fraturas do úmero. A conexão rígida entre os parafusos de bloqueio e a placa permite melhores resultados clínicos em fraturas simples, fraturas cominutas e em ossos osteoporóticos da região proximal do úmero.



Its indication, purpose and/or use is the stabilization and fixation of different types of fractures, providing the recovery of humeral fractures. The rigid connection between the locking screws and the plate allows for better clinical results in simple fractures, comminuted fractures and in osteoporotic bones in the proximal region of the humerus.



Su indicación, finalidad y/o uso es la estabilización y fijación de diferentes tipos de fracturas, proporcionando la recuperación de las fracturas humerales. La conexión rígida entre los tornillos de bloqueo y la placa permite mejores resultados clínicos en las fracturas simples, fracturas conminutas y en huesos osteoporóticos en la región proximal.



ANVISA 80062900021

Aço Inoxidável (NBR ISO 5832-1)

Stainless Steel (NBR ISO 5832-1)

Acero Inoxidable (NBR ISO 5832-1)



Biologicamente compatível,
evitando reação inflamatória.

Biocompatible, avoiding
inflammatory reactions.

Biologicamente compatible,
evitando reacción inflamatoria.



Formato anatômico

Anatomic shape

Formato anatómico



Acabamento polido

Polished finish

Acabado pulido



Furos de bloqueio: $\varnothing 4,0\text{mm}$,
e simples: $\varnothing 4,5\text{mm}$

Locking holes: $\varnothing 4,0\text{mm}$,
and simple: $\varnothing 4,5\text{mm}$

Perforaciones de bloqueo:
 $\varnothing 4,0\text{mm}$ **y simples:** $\varnothing 4,5\text{mm}$



Utilizado com parafusos:

Bloqueio: $\varnothing 3,5\text{mm}$
Cortical: $\varnothing 3,5\text{mm}$
Esponjoso: $\varnothing 4,0\text{mm}$

Use with screws:

Locking: $\varnothing 3,5\text{mm}$
Cortical: $\varnothing 3,5\text{mm}$
Cancellous: $\varnothing 4,0\text{mm}$

Utilice con tornillos:

De bloqueo: $\varnothing 3,5\text{mm}$
Cortical: $\varnothing 3,5\text{mm}$
Esponjoso: $\varnothing 4,0\text{mm}$

PLACA BLOQUEADA LPS ÚMERO PROXIMA EM AÇO

LOCKING PLATE LPS PROXIMAL
HUMERUS IN STAINLESS STEEL

PLACA BLOQUEADA LPS HÚMERO
PROXIMAL EN ACERO



ANVISA 80062900021

CÓDIGO <i>CODE</i> CÓDIGO	Nº FUROS <i>Nº HOLES</i> Nº PERF.	LARGURA (MM) <i>WIDTH (MM)</i> ANCHO (MM)	COMPRIM. (MM) <i>LENGHT (MM)</i> LONGITUD (MM)	ESPESSURA (MM) <i>THICKNESS (MM)</i> ESPEJOR (MM)
5711-004-003	3	12	95	3,0
5711-004-004	4	12	110	3,0
5711-004-005	5	12	125	3,0
5711-004-006	6	12	140	3,0
5711-004-007	7	12	155	3,0
5711-004-008	8	12	170	3,0
5711-004-009	9	12	185	3,0
5711-004-010	10	12	200	3,0
5711-004-011	11	12	215	3,0
5711-004-012	12	12	230	3,0